

INSTRUCTIONS FOR USE

IFU CODE : 03



E60 – E63 Sach Foot



E61 Single Axis Foot



E62 Dynamic Foot



INFORMATION



Last Update: 01.06.2026

- Please read this document carefully before using the product.
- Instruct the user in the appropriate and safe use of the product.
- Follow the safety instructions to abstain injuries and damage to the product.
- Please keep this document in a safe place.

Contents

1. E60 – E63 Sach Foot	1
2. E62 Dynamic Foot	2
3. E61 Single Axis Foot and Adapters	3
4. Environmental Conditions	4
5. General Safety Instruction	5
6. Cleaning	5
7. Warranty	6

1. E60 – E63 Sach Foot

Sach foot has split toes.

Enclosed heel bumper.

Two different color options available (Beige and Brown).



Sach foot is used with E01 (foot adapter) and E60 mika. Torque wrench should be used to tighten screws. In order to set the wrench, following power should be applied.

Torque [Nm]	15 Nm (E01s-M8, E01T-M8)
	30 Nm (E01s-M10, E01T-M10)

1.1 Properties

Foot Size	22-29
Load Weight	100 kg
Heel Height	10 mm
Color	Beige-Brown
Side	L/R

2. E62 Dynamic Foot

Manufactured with soft heel material and belt inside.

Dynamic foot has split toes.

Dynamic foot is used with E01 (foot adapter) and E 60 mika.

Torque wrench should be used

to tighten screws. In order to set the wrench, following power should be applied.



Torque [Nm]	15 Nm (E01s-M8 & E01T-M8)
	30 Nm (E01s-M10 & E01T-M10)

2.1 Properties

Foot Size	22-29
Load Weight	100 kg
Heel Height	10 mm
Side	L/R



NOTICE

This product is waterproof. After exposure to chlorinated or saltwater, it is critical to rinse the device thoroughly with clean fresh water, drain the water by holding it vertically, and dry it with a lint-free cloth to maintain its mechanical lifespan.

3. E61 Single Axis Foot

The Single Axis Foot is suitable for lower extremity prosthesis applications.

The Single Axis Foot has smooth surface and natural shape.



Single axis foot is used with E50 (Modular single axis foot adapter) and E61 Cup. Torque wrench should be used to tighten screws. In order to set the wrench, following torque should be applied.

Torque [Nm]	6 Nm (E50)
--------------------	------------

3.1 Properties

Foot Size	22-29
Maximum Body Weight	100 kg
Heel Height	10 mm
Side	L/R

4. Environmental Conditions

E62 Dynamic Foot

Acceptable Environmental Conditions

Fresh water, soapy water, chlorinated water, sea water

Temperature range: -10 °C to +45°C

Moisture: Allowable relative humidity 20 % to 90 %, non-condensing

Unsuitable Environmental Conditions

Chemicals/liquids: Urine, acids e.g.

Solids: Dust in high concentrations, highly hygroscopic particles (e.g., talcum), intensive contact with sand

No impacts or mechanical vibrations

E60 – E63 Sach Foot, E61 Single Axis Foot

Acceptable Environmental Conditions

Temperature range: -10 °C to +45°C

Moisture: Allowable relative humidity 20 % to 90 %, non-condensing

Unsuitable Environmental Conditions

Chemicals/liquids: Water, perspiration, urine, acids e.g.

Solids: Dust in high concentrations, highly hygroscopic particles (e.g., talcum), intensive contact with sand

No impacts or mechanical vibrations

5. General Safety Instruction



CAUTION!

To avoid the risk of injury and product damage, do not use for longer than the service life of the product.

The prosthetic foot should not be sanded.

Keep the product away from flammable materials. Protect the product against heat.

Avoid using this product in unsuitable environmental conditions. If the product has been exposed to unsuitable environmental conditions, check for damage.

Do not use the product that is damaged, has limited functionality or is in doubtful condition. Ensuring that appropriate measures are taken (e.g. control, cleaning, repair, replacement by the authorized workshop or manufacturer).

To prevent mechanical damage, be careful and caution when working the product.

Use the product for a single patient to avoid injuries and product damage.

Products should be used in appropriate combinations. Using the product with unsuitable combinations may result in possible injury.

6. Cleaning

- 1.Clean the product with a moisture, soft cloth.
- 2.Then dry the product with a soft cloth
- 3.Afterwards, the remaining moisture should be left to air-dry.

7. Warranty

The IB-ER company provides warranty for the products in case of material and production flaws. Warranty duration for the above-mentioned products is individually mentioned below.

Parts	Duration after the IB-ER Invoice Date
E60	6 Months
E61	
E62	
E63	

This warranty only covers the replacement or repair of failed products, provided that assembly was performed by qualified individuals or authorized institutions. This warranty does not apply under the following conditions:

- Improper use of the product, use under unsuitable conditions, or faulty assembly.
- Issues arising from the replacement of our products with third-party components.



İB-ER Makina A.Ş
Selvilitepe OSB Mah. OSB 2006. Cad. No: 13/1
45400, Turgutlu - Manisa / Turkey

Phone : +90 (236) 312 32 00
Fax : +90 (236) 312 32 06
info@ib-er.com
www.ib-er.com



Ercan Gozuacik
Reinickendorfer Str 45 13347 Berlin GERMANY
(for regulatory affairs only)